

## IE-11 (defender IPL/SHR) Parameters setting table

### Note:

The IPL energy in above table is just the starting value, if you feel nothing with this energy, please raise it one by one gradually.

Because every people have different skin color, skin bear capability, so when treatment, the rules is from low energy to high energy, always ask customers feeling, to find the best suitable energy (the customer can feel hot energy within their bear capability)



### 1.SHR

Please shave hair 1-2 days in advance, Increase energy based on customer experience, if smells the scorched hair, removes the hair with tweezers, no pain, the best effect. One treatment per month, must be continued 3-4 times, and the course of treatment will be increased until the permanent epilation result is achieved

Skin Type	Pulse Width	Frequency	Energy(J)	Cooling
I	4-5	3-5	14~28	2-3
II	4-5	3-5	12~25	2-3
III	3-4	3-4	10~23	2-3
IV	3-4	3-4	8 ~20	3-4
V	2-3	2-3	6~15	3-4



- 1, choose skin color
- 2, Choose parts of body to treat
- 3, Hair Status
- 4, Adjust Energy “10-28J/cm<sup>2</sup>” “<” mean : Energy reduce ; Energy “>” mean :Energy increase
- 5, Adjust :4-5 pieces”therapeutic temperature of the handle”
- 6, Adjust Frequency, usually“3-5HZ”, “<” reduce, “>”increase

After adjust parameters(skin color, treatment part, energy, Frequency etc. Please press  to  , keep about 3-5 second, then parameters will be stored, you will not need adjust when you use it next time.

Friendly reminder:(Parameters need to be adjusted appropriately according to treatment part of body and customer experience)

- 1)Energy :13-28J/cm<sup>2</sup>
- 2)Snowflake: 3-3.5/piece
- 3)Frequency: 3-5HZ (Light output speed)
- 4)Pulse width: 5-8Ms
- 5)Skin color:I, II, III, IV, V.
- 6)Parts of body : Face, arm, body, Bikini, Leg