


Multi11 settings table




After set settings ready, please press  to save settings. Otherwise, the settings will not be saved.

1. SHR

Skin Type	Pulse Width	Frequency	Energy	Cooling
White Skin	7	3-5	From 15J	4-5
Yellow Skin	4-5	3-4	From 15J	4-5
Dark Skin	Not recommended to do !!!			

2. IPL



A: Hair Removal (640nm)

Skin Type	Pulse Width	Pulse Interval	Pulse Number	Energy	Cooling
White Skin	4-5	5-6	5-6	From 10J	4-5
Yellow Skin	3-4	6-7	4-6	From 7J	4-5
Dark Skin	Not recommended to do !!!				

B: Spider Veins Removal (590nm)

Skin Type	Pulse Width	Pulse Interval	Pulse Number	Energy	Cooling
White Skin	4-5	5-6	2-3	From 10J	4-5
Yellow Skin	3-4	6-7	2-3	From 5J	4-5
Dark Skin	Not recommended to do !!!				

C: Skin Rejuvenation (530nm)

Skin Type	Pulse Width	Pulse Interval	Pulse Number	Energy	Cooling
White Skin	4-5	5-6	3-4	From 10J	4-5
Yellow Skin	3-4	6-7	3-4	From 5J	4-5
Dark Skin	Not recommended to do !!!				

D: Speckle Removal (530nm)

Skin Type	Pulse Width	Pulse Interval	Pulse Number	Energy	Cooling
White Skin	4-5	5-6	3-5	From 10J	4-5
Yellow Skin	3-4	6-7	3-5	From 5J	4-5
Dark Skin	Not recommended to do !!!				

E: Acne Removal (430nm)

Skin Type	Pulse Width	Pulse Interval	Pulse Number	Energy	Cooling
White Skin	4-5	5-6	2-3	From 6J	4-5
Yellow Skin	3-4	6-7	2-3	From 3J	4-5
Dark Skin	Not recommended to do !!!				

3. YAG

Item	Energy	Frequency	Distance
Eyebrow Removal	500-600mJ	2~3HZ	1~2cm

Tattoo Removal	700-800mJ	3~5HZ	3~5cm
Nevus Removal	700-800mJ	3~5HZ	4~5cm
Birthmark Removal	1000-1200mJ	3~5HZ	4~5cm
Nevus of Ota	1000-1200mJ	3~5HZ	4~5cm
Skin Rejuvenation	700-800mJ	3~5HZ	4~5cm

Note:

The IPL energy in above table is just the starting value, if you feel nothing with this energy, please raise it one by one gradually.

Because every people have different skin color, skin bear capability, so when treatment, the rules is from low energy to high energy, always ask customers feeling, to find the best suitable energy (the customer can feel hot energy within their bear capability)

Regarding the results, please be attention to below:

Step 1: how to find good settings. Different people needs different settings sometimes as everyone has different skin conditions or skin colors.

The principle to find good settings is:

try a low settings firstly on customer's skin(you can use the settings we recommended in the table) and ask about feeling. If she feels nothing or feels very weak, increase energy, until she feels "more than warm but not very heat".

Step 2: after finding good settings in Step 1, use this good settings to treat, until "skin become slight red". Then the treatment this time ends.

How fast I slide handle on skin? It depends on what frequency you set on skin.

The next shot should be cover a little bit on the previous shot, in order to avoid any skin area where is not fired at.

According to our customers' feedback to produce the best results and keep the machine's best performance, not increase the frequency more than 5Hz.

Remark: some customer may feel sensitive to this light energy even the energy is lower than other normal customers. In this situation, you can choose a way: low energy but more firing times in one treatment, until skin become slightly red.

Attention: skin must receive energy gradually and results can be good. Never use too high energy in order to make skin slightly red to end treatment.

The method above also applies to IPL treatments.

And when to use SHR handle sliding on skin, don't slide too fast, or the energy cannot cover treatment area evenly and it will cause that some areas on skin absorb more energy but some others absorb less or nothing, and finally hair on some area may stop growing but still grow on other area.